

LABCOLD SPARKFREE CHEST FREEZER- 607L

125-052

**PRODUCT DESCRIPTION**

The Labcold Sparkfree Chest Freezer provides improved temperature uniformity and stability with much-reduced temperature rise following door openings

Labcold Sparkfree Freezers are designed specifically for the safe storage of flammable materials in a laboratory environment. All sources of ignition, such as the thermostat are removed from the inside of the cabinet.

Key Features:

- Sparkfree Interior
- Door Lock
- Easy-Grip Door Handle
- CFC & HCFC Free
- Temperature Set-Point Range: - 18 to - 22 C
- Visual High-Temperature Alarm
- Operator Controlled Defrost

Specifications:

| | |
|------------------|---------------|
| Capacity: | 607 Litres |
| Style: | Chest Freezer |

| | |
|---------------------------------|---------------------|
| Exterior Dimensions(mm): | 860H x 1700W x 730D |
| Interior Dimensions(mm): | 720H x 1580W x 530D |
| Weight: | |
| Energy Consumption: | 1.01kWh/24hr |
| Operating Temperature: | -18 C to - 25 C |
| Factory Set Point: | -20 C |
| Interior Construction: | Aluminium |
| Exterior Construction: | White Enamel |
| Lid Style: | Solid |
| Baskets | 2 |
| | |
| Feet: | 4 (adjustable) |
| High Temperature Alarm: | Visual |

ADDITIONAL INFORMATION

| | |
|------------------------|---------|
| Animal | N/A |
| SKU | 125-052 |
| Colour/Style | White |
| Brand | Labcold |
| Alarm | Yes |
| Fridge Capacity | 607L |
| Door Type | Regular |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights