

**LABCOLD SPARKFREE CHEST FREEZER- 215L**

125-050

**PRODUCT DESCRIPTION**

**The Labcold Sparkfree Chest Freezer provides improved temperature uniformity and stability with much-reduced temperature rise following door openings**

Labcold Sparkfree Freezers are designed specifically for the safe storage of flammable materials in a laboratory environment. All sources of ignition, such as the thermostat are removed from the inside of the cabinet.

**Key Features:**

- Sparkfree Interior
- Door Lock
- Easy-Grip Door Handle
- CFC & HCFC Free
- Temperature Set-Point Range: - 18 to - 22 C
- Visual High-Temperature Alarm
- Operator Controlled Defrost

**Specifications:**

<b>Capacity:</b>	215 Litres
<b>Style:</b>	Chest Freezer

<b>Exterior Dimensions(mm):</b>	860H x 725W x 730D
<b>Interior Dimensions(mm):</b>	720H x 580W x 530D
<b>Weight:</b>	46
<b>Energy Consumption:</b>	1.01kWh/24hr
<b>Operating Temperature:</b>	-18 C to - 25 C
<b>Factory Set Point:</b>	-20 C
<b>Interior Construction:</b>	Aluminium
<b>Exterior Construction:</b>	White Enamel
<b>Lid Style:</b>	Solid
<b>Baskets</b>	1
<b>Basket Dimensions</b>	200H x 250W x 500D
<b>Feet:</b>	4 (adjustable)
<b>High Temperature Alarm:</b>	Visual

## ADDITIONAL INFORMATION

<b>SKU</b>	125-050
<b>Overall Product Dimensions (mm)</b>	860H x 725W x 730D
<b>Colour/Style</b>	White
<b>Brand</b>	Labcold
<b>Alarm</b>	Yes
<b>Fridge Capacity</b>	215L
<b>Door Type</b>	Regular

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights