

**CONTROLLED DRUGS CABINET - 335W X 270D X 300HMM**

115-491

**PRODUCT DESCRIPTION****Single point locking controlled drugs cabinets that exceed the statutory requirement for the storage of controlled drugs**

Ensuring security with a strong 1.5mm (16 gauge) mild steel, seamless welded construction and fitted with a BS3621 anti-pick & anti-drill, six-pin cylinder key retaining lock. Supplied with 2 working keys, single shelf and identity card.

**Key Features:**

- Supplied right hand hinging as standard, left hand available on request
- Single-point lock
- Supplied with 2 Keys & ID Card
- 1 x Shelf
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs per profile
- Lighting transformer casing adds 105mm to the overall height
- 1.5mm (16 gauge) mild steel, seamless welded construction
- Mounting holes in back (fixings not included)
- Adjustable & removable shelves
- Modular design allows for easy interchange with HTM 63 storage systems
- Painted white finish
- Compatible with ABLOY® CLIQ™ Remote, integrating electronics into high-security mechanical disc cylinders & keys, combining electronic security & remote management in real time. Can be supplied with a micro-switch to allow connection to remote lighting/nurse call system by others (CD905), or supplied with a 12v DC LED lighting/warning system with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC

### Dimensions:

- 335W x 270D x 300Hmm
- 9.0kg

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	115-491
<b>Overall Product Dimensions (mm)</b>	335W x 270D x 300Hmm
<b>Colour/Style</b>	White
<b>Brand</b>	Bristol Maid

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights