

CONTROLLED DRUGS CABINET - 500W X 450D X 1250H MM

115-496



PRODUCT DESCRIPTION

Controlled drugs cabinet that exceed the statutory requirement for the storage of controlled drugs

Ensuring security with a strong 1.5mm (16 gauge) mild steel, seamless welded construction and fitted with twin BS3621 anti-pick & anti drill, six-pin cylinder key retaining locks. Supplied with three shelves and identity card.

Key Features:

- Supplied right hand hinging as standard, left hand available on request
- Single lock
- Supplied with 2 Keys & ID Card
- 4 x Shelves
- Key profile unique to Hospital Metalcraft Ltd, offering 250,000 differs per profile
- Lighting transformer casing adds 105mm to the overall height
- 1.5mm (16 gauge) mild steel, seamless welded construction
- Mounting holes in back (fixings not included)
- Adjustable & removable shelves
- Modular design allows for easy interchange with HTM 63 storage systems
- Painted white finish
- Compatible with ABLOY® CLIQ™ Remote, integrating electronics into high-security mechanical disc cylinders & keys, combining electronic security & remote management in real time. Can be supplied with a micro-switch to allow connection to remote lighting/nurse call system by others (CD905), or supplied with a 12v DC LED lighting/warning system with volt-free contact to allow connection to light/nurse call systems. Requires 100v - 240v ~, 50Hz -60Hz, 0.3A Mains supply. Supplied fitted with a UK 3 pin plug. CE marked & conforms to the low voltage directives 73/23/EEC & 89/336/EC

Dimensions:

- 500W x 450D x 1250H mm
- 59.0kg

ADDITIONAL INFORMATION

Animal	N/A
SKU	115-496
Overall Product Dimensions (mm)	500W x 450D x 1250H mm
Colour/Style	White
Brand	Bristol Maid

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights