

## HANDS FREE WASTE BIN WHITE LID (FOR GENERAL USE)

225-45-CONFIG-WHITE



### PRODUCT DESCRIPTION

## Flame retardant construction handfree bin, compliant with HTM07-01 (Safe Management of Healthcare Waste) standard set by the NHS

White Lid version for General waste, Waste for incineration (Yellow lid) available separately.

#### Product Features

- Hands free operation, foot pedal must be used to open
- Silent - dampened closing lid as standard
- Sealed base prevents leakages & facilitates easy cleaning
- Integrated wall mounting bracket as standard
- Designed and tested to HTM 83
- Sizes HxWxD: 20 litre = 51 x 32 x 33cm, 50 litre = 59 x 42 x 43cm, 70 litre 78 x 42 x 43cm

### ADDITIONAL INFORMATION

|  |                     |
|--|---------------------|
| <b>SKU</b>                             | 225-45-config-white |
| <b>Brand</b>                           | Sunflower Medical   |
| <b>Animal</b>                          | N/A                 |
| <b>Overall Product Dimensions (mm)</b> | Various             |
| <b>Colour/Style</b>                    | White               |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights