

## HANDS FREE WASTE BIN YELLOW LID (FOR INCINERATION)

225-45-CONFIG-YELLOW



### PRODUCT DESCRIPTION

## Flame retardant construction handfree bin, compliant with HTM07-01 (Safe Management of Healthcare Waste) standard set by the NHS

Yellow Lid version for Waste for incineration, recycling (white lid), and domestic waste (black lid) available separately.

#### Product Features

- Hands free operation, foot pedal must be used to open
- Silent - dampened closing lid as standard
- Sealed base prevents leakages & facilitates easy cleaning
- Integrated wall mounting bracket as standard
- Designed and tested to HTM 83
- Sizes HxWxD: 20 litre = 51 x 32 x 33cm, 50 litre = 59 x 42 x 43cm, 70 litre 78 x 42 x 43cm

### ADDITIONAL INFORMATION

<b>SKU</b>	225-45-config-yellow
<b>Brand</b>	Sunflower Medical
<b>Animal</b>	N/A
<b>Overall Product Dimensions (mm)</b>	Various
<b>Colour/Style</b>	Yellow

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights