

## SURGEONS DECK MOUNTED 6" LEVER TAP

115-250-3



### PRODUCT DESCRIPTION

## Ideal for use on surgical sinks

Finished in a high shine chrome, this high quality tap is fitted with 6" levers and ceramic ¼ turn seats for rapid on/off flow

#### Features

- Swivelling spout and divided flow
- Offset mountings enable fitting on any sink / surface with holes centres between 160mm and 240mm.
- Inlet size 1/2" BSP
- Deck Mounted

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	115-250-3
<b>Colour/Style</b>	Stainless Steel
<b>Operation</b>	Surgeon (Lever)

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights