

A4 CAGE CLIP BOARD

165-400-1290

**PRODUCT DESCRIPTION**

A4 Portrait Stainless steel clipboard that is designed to last a life-time, securely holding patient information to your cage doors

- Made from a durable & hygienic T504 stainless steel allowing the clipboards to be easily wiped down and disinfected.
- Product includes bracket to fit a standard mesh cage and kennel door.
- Designed for Burtons cages

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	165-400-1290
Colour/Style	Stainless Steel
Brand	Burtons
Fit	Cage
Type	Bank Accessories

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights