

## SCISPIN MINI CENTRIFUGE

146-032-COLOURS



### PRODUCT DESCRIPTION

---

**The SciSpin mini centrifuge is great for practices with limited space**

#### Key Features:

- Ultra compact design
- Max 7,000 RPM /2,680g
- Includes 2 rotors - the first spins at 8 x 1.5 or 2.2ml blood tubes and the second 16 x 0.2ml PCR tubes (used for DNA testing)
- Dual lid lock for safety
- Ultra quiet <45dB
- Perfect for quick spin downs and microfiltration
- Rapid acceleration and deceleration
- Very easy to use with no set programmes - simply stop and start by opening and closing the lid
- Dimensions: 150mmW x 150mmD x 117mmH
- Weight: Just 0.5kg

#### Specifications

- Rotor Capacity: 8 x 1.5/2.2ml
- 16 x 0.2ml PCR
- Max Speed: 7,000 RPM
- Max RCF: 2, 680g

- Timer: Continuous Operation
- Noise Level: 45db
- Power 100-240V 50/60Hz

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	146-032-colours
<b>Generic Size</b>	X-Small
<b>Catagory Sub Heading</b>	Colour Options Available

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights