

SWANN MORTON STERILE STITCH CUTTER

400-073



PRODUCT DESCRIPTION

Traditional carbon steel sterile standard stitch cutter

With its rounded bull nose tip which ensures there is no damage to the skin whilst the surgical edge cuts through the suture with minimal resistance.

Convenient and cost effective these sterile stitch cutters are ideal for use in 'Out Patient clinics', within 'Accident and Emergency departments', on ' Hospital Wards' in 'General Practice' and for 'Nursing in the Community'.

- Pack of 100

ADDITIONAL INFORMATION

Animal	Mixed
SKU	400-073
Brand	Swann Morton

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights