

TERUMO 3 PART DISPOSABLE SYRINGE WITH NEEDLE

400-01-WITHNEEDLE



PRODUCT DESCRIPTION

Quality Disposable syringes for precise administration of drugs

Key Features

- For the injection of fluids
- Can also be used for aspiration of fluids, the sampling of blood, or other body fluids
- Latex free gasket
- Bold markings for ease of measurements and dosing
- Individually sterile wrapped
- Constructed of polypropylene, polystyrene, thermoplastic elastomer & silicone oil

ADDITIONAL INFORMATION

Animal

Mixed

SKU

400-01-withneedle

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights