

## TERUMO 3 PART DISPOSABLE SYRINGE

400-00-WITHOUTNEEDLE



### PRODUCT DESCRIPTION

---

## Quality Disposable syringes for precise administration of drugs

#### Key Features

- For the injection of fluids
- Can also be used for aspiration of fluids, the sampling of blood, or other body fluids
- Latex-free gasket
- Bold markings for ease of measurements and dosing
- Individually sterile wrapped
- Constructed of polypropylene, polystyrene, thermoplastic elastomer & silicone oil

#### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	400-00-withoutneedle

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights