

## BURTONS CAT CUBE BOARDING KENNEL

180-026-BUNDLE



### PRODUCT DESCRIPTION

---

**The new construction of The Cat Cube uses our newest material that is anti-bacterial, and completely water resistant. With toughened glass doors as standard. Giving your feline patients a warm, safe, quiet & comfortable place to stay**

Cat Cubes, also known as Condos, have been designed around feline needs, with features such as a built-in resting shelf and pass through horizontal portholes that can be opened or closed allowing increased space for the patient or access to litter.

Our new Cat Cubes are made with a foam core of closed cell PVCU with a solid PVC surface. This specifically chosen board is non-toxic, formaldehydes free, non-porous and resistant to bacteria and completely impervious to liquid. Which means you get a cat cube that is comfortable, easy to clean and long lasting! Not to mention that WCB board is also made completely from recycled materials and can also be recycled too!

Cat cubes use a simple stack design, making installation easy and can be used in conjunction with a purpose made laminate cat pedestal raising the bank to a comfortable working height and giving extra storage in the boarding area.

The pass through horizontal portholes work in conjunction with external sliding dividers which can be opened or closed from outside the cube.

#### Doors

The toughened glass doors allow for easy patient observation with ample ventilation and a simple one-touch push, quiet catch.

## Banks

Cat cubes are sold in left, middle and right hand designs to allow for correct pass through set up, if you require a smaller 2 cube bank simply miss out the middle cube, or if you require a longer run just add more middle cubes to your order. Cat cubes require a minimum size of 2x1 bank with the ideal bank being 3x2

All of our cages are available in standard sizes from stock with special made to measure sizes available on request. The completely modular, standardised design ensures your bank can grow with your practice.

If hospitalisation is required, we would recommend using our 'Quiet Time' Polypropylene or 'Lifetime' Stainless Steel veterinary cages, which are specifically designed for use in a clinical care environment.

## Features

- Cosy & quiet atmosphere
- Vented, toughened glass doors for advanced patient observation
- Horizontal portal with external sliding dividers
- Laminate plinth and power bar options
- Robust, anti-bacterial & waterproof construction
- Built-in resting shelf
- Bespoke sizes available
- One-touch latch lock
- Modular & stackable design
- Easy to install and set up
- Minimum 1x2 design but no maximum
- Standard Size: 24" W x 30" H x 26.4 D

## Let us design it for you...

At Burtons, we offer a completely complimentary custom design service, meaning whatever your plan for the room, whether it is to gain ISFM certification or just maximising every bit of space. Our dedicated projects team will use their expert knowledge to complete a design that will be both functional and to specification.

We can arrange for a full site survey free of charge or you can simply provide us with your room dimensions and our projects team will complete a free 3D room design based around your needs.

## Installation

If it's a single cage or a new practice setup the easy installation process is straight forward and fast. Once your happy with your order, all our cages are delivered pre-made and ready to go, but if you prefer why not take advantage of our in-house installation team, who can then arrange a suitable time to build your new bank, working with you to minimize disruption.

## Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been

engineered from the ground up to last!

### Care of your Burtons Cat Cubes

- Applying cleaner directly to the surface of the cube can cause staining— apply using cleaning cloth instead.
- Always dilute the cleaner to the manufacturer’s guidelines.
- Dry up any standing liquid
- Abrasive cleaning materials should not be used.

### Maintenance of your Laminate Burtons Cat Cubes

From time to time dependent upon cleaning regimes and animal behaviour the silicon sealant may need to be renewed, this is **not** covered by Burtons unlimited first-year warranty as it is considered to be a consumable item. It is recommended that a high modulus silicone is used to reseal the units. High modulus sealants withstand cleaning much better and make a stiffer and more rigid seal. A competent handyman can renew the sealant. This will not affect your manufacturer’s warranty.

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	180-026-bundle
<b>Overall Product Dimensions (mm)</b>	24in x 30in x 26.5in
<b>Colour/Style</b>	White
<b>Brand</b>	Burtons
<b>Catagory Sub Heading</b>	New design!
<b>Type</b>	Kennel
<b>Depth</b>	26.5"
<b>Complete Bank</b>	No
<b>Material</b>	WCB Board
<b>Mobile</b>	No
<b>Door Type</b>	Toughened Glass

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights