

OXYGEN THERAPY DOOR FOR BURTONS CAGE

165-450-BUNDLE



PRODUCT DESCRIPTION

Transform your cage into an oxygen therapy chamber

High quality oxygen therapy door featuring a thermometer/hygrometer mounted on the door, which monitors cage temperature and humidity for optimal treatment.

The door fits onto the front of a Lifetime or Quiet time cage in place of the regular door, A rubber gasket seal is open at the corners to uniformly distribute oxygen and exhaust carbon dioxide.

Features

- Uses regular hinges and latches to close securely
- Additional port is provided to allow the fitting of an oxygen analyser
- Comes with humidifier and ice tray.

ADDITIONAL INFORMATION

| | |
|-----------------------------|-------------------------|
| Animal | Small Animal |
| SKU | 165-450-bundle |
| Colour/Style | Stainless Steel |
| Brand | Burtons |
| Catagory Sub Heading | Various Sizes Available |
| Fit | Cage |
| Type | Bank Accessories, Cage |

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights