


## ACRYLIC DOOR FOR BURTONS CAGE

165-480-\*BUNDLE



### PRODUCT DESCRIPTION

**Taking the place of your standard door the acrylic door allows for much easier patient monitoring and helps reduce the spread of disease**

#### Features

- Constructed from 6mm clear acrylic
- Built in vent holes
- Stainless steel hinge and latch mechanism
- Simple to fit with a full range of sizes available

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	165-480-*bundle
<b>Colour/Style</b>	Clear
<b>Brand</b>	Burtons
<b>Catagory Sub Heading</b>	Various Size Avaliable
<b>Type</b>	Bank Accessories, Cage

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights