

CAGE HINGE IV HOOK FOR BURTONS CAGES

165-410-001-BUNDLE



PRODUCT DESCRIPTION

Stainless steel hinge pin IV hooks replace the cage hinge pin to provide convenient and secure positioning of an IV bag

- Removes potential trip hazards from kennel area
- Easy to attach
- Fits any Burtons cage

ADDITIONAL INFORMATION

| | |
|---------------------|--------------------|
| Animal | Mixed |
| SKU | 165-410-001-bundle |
| Colour/Style | Silver |
| Brand | Burtons |
| Fit | Cage |
| Type | Bank Accessories |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights