

## CAGE HINGE IV HOOK FOR BURTONS CAGES

165-410-001-BUNDLE



### PRODUCT DESCRIPTION

**Stainless steel hinge pin IV hooks replace the cage hinge pin to provide convenient and secure positioning of an IV bag**

- Removes potential trip hazards from kennel area
- Easy to attach
- Fits any Burtons cage

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	165-410-001-bundle
<b>Colour/Style</b>	Silver
<b>Brand</b>	Burtons
<b>Fit</b>	Cage
<b>Type</b>	Bank Accessories

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights