

## JONES CROSS ACTION TOWEL

400-107



### PRODUCT DESCRIPTION

## Jones Cross Action Towel Clamps 9cm (3.5in)

Designed for securely fixing drapes during procedures

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	400-107
<b>Left/Right Handed</b>	Right
<b>Tissue Forcep Type</b>	Jones Cross Action
<b>Instrument Size</b>	9cm

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

