

## HARVEY ELITE STETHOSCOPE (28IN)

WA5079-CONFIG



### PRODUCT DESCRIPTION

---

## The Harvey Elite stethoscope is equipped with dual-bore, latex-free tubing that creates two uninterrupted sound channels

The Elite's diaphragm is designed to accentuate important sounds all this enables you to hear cardiac, pulmonary and vascular sounds.

- Full-range cardiology / pulmonology stethoscope
- Double-head (bell and flat diaphragm)
- Stainless steel chestpiece with nonchill rim
- Rotatable stainless steel binaurals and interchangeable comfort sealing eartips
- Triple-leaf binaural spring encased in polyurethane allows size adjustment while preventing breakage
- 28" length- Latex-free
- Choice of 4 Colours

### ADDITIONAL INFORMATION

**Animal**

Small Animal

**SKU**

WA5079-config

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights