

**DARAY MAGW150 HANDHELD WOODS LIGHT**

145-856

**PRODUCT DESCRIPTION**

**The MAGW150 is specifically designed to meet the demanding needs of today's veterinary practices whilst providing the finest quality design with superior performance, reliability and value**

- The MAGW150 emits ultraviolet light with the ability to note minor skin colour changes and variations in skin tone.
- Under UV light different substances become luminous when exposed making it a useful tool for veterinary surgeons.

**ADDITIONAL INFORMATION**

|               |         |
|---------------|---------|
| <b>Animal</b> | N/A     |
| <b>SKU</b>    | 145-856 |

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

