

## MACINTOSH LARYNGOSCOPE BLADES

117-6-MAC-BUNDLE



### PRODUCT DESCRIPTION

## High quality stainless steel Macintosh blades

The Macintosh Laryngoscope Blade is the predominate model of the curved blade types. During laryngoscopy, the tip or beak of the blade is compressed into the angle formed by the base of the tongue and the epiglottis, indirectly raising the epiglottis. Some sizes are available flangeless

#### Key Features:

- Stainless steel construction and a satin non-glare finish
- Blades compatible with all major brands of standard illumination laryngoscopes
- Multiple sizes available

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	117-6-mac-bundle
<b>Colour/Style</b>	Stainless Steel
<b>Catagory Sub Heading</b>	Range of Sizes Available

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights