

## MILLER LARYNGOSCOPE BLADES

117-6-MIL-BUNDL



### PRODUCT DESCRIPTION

## High-quality Miller Laryngoscopes blades

The Miller Laryngoscope Blade is the most popular of the straight laryngoscopes. The side of the flange is reduced to minimize trauma and the curve at the beak or tip is extended to improve lifting of the epiglottis. These improvements facilitate greater exposure of the larynx in difficult-to-intubate patients.

- Blade flattened to move the patient's tongue aside, rather than compressed
- Flat base to help to hold the tongue out of the field of view
- Overall field of view increased by modifying the cross-sectional area of the blade without altering width or height

#### Key Features

- Stainless steel construction
- Blades compatible with all major brands of standard illumination laryngoscopes
- Multiple sizes available

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	117-6-mil-bundl
<b>Colour/Style</b>	Stainless Steel
<b>Catagory Sub Heading</b>	Range of Sizes Available

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights