

## BURTONS RAIL MOUNTED ANAESTHESIA MACHINE

102-088-BUNDLE



### PRODUCT DESCRIPTION

---

## Fully utilise medical rail systems with this specially designed Burton's anaesthetic machine!

Ideal for piped gas systems this medical rail mounted, O<sub>2</sub> only, anaesthetic machine is a perfect addition to any practice. The machine is built to last, with a stainless steel frame and Burton's own high-quality components. Supplied as standard with a high accuracy rotameter and emergency O<sub>2</sub> flush (including pressure release safety valve).

Installation and moving the system is easy on your medical rail system, just clip on and go, with built-in lock down points to stop any movement whilst the system is in use.

Use this machine alongside a piped gas system to instantly save precious floor space and remove the need to move large cylinders around the practice.

#### All models are supplied with :

- 1-metre gas supply lead with indexed BS Schrader probes
- High accuracy rotameters for precise control of gas flow rates
- Corner unit outlet with built-in emergency oxygen flush
- Non-return valve and safety valve
- Stainless steel frame

#### Don't forget your Vaporiser

The Series 3 is a temperature compensating, out of circuit vaporiser in a traditional design, familiar to most veterinary practitioners. Calibrated specifically for the specified anaesthetic drug using a laser refractometer.

#### Fitting Type

Choose from the standard taper mounted fitting (cage mount) (UNI) or the quick release (QR) option; allowing you to easily swap the vaporiser to different anaesthetic machines throughout the practice. Additional quick release compatibility blocks will be required for each anaesthesia machine. (Available separately)

### Key filled or Funnel (Pour fill)

Key filled (KF) systems require the addition of key filler adaptors but are generally safer and have a lot less wastage compared to standard funnel filled (FF) vaporisers, saving you money in the long run.

### Technical Specification:

- O<sub>2</sub> only
- Overall size: 21cmD x 40cmW x 47cmH
- Outlet connector: BS22mm Male
- Flow control: 0.1 to 10 l/min

### Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

*Please Note: Vaporiser Sold Separately*

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	102-088-bundle
<b>Colour/Style</b>	Stainless Steel
<b>Gas Type</b>	O <sup>2</sup>
<b>Gas Supply</b>	'E' Size Cylinder (s), 'F' Size Cylinder(s), Piped Gas
<b>Flowmeter Type</b>	Graduated Rotameter
<b>Concentrator</b>	No
<b>Optional O2 Alarm</b>	Yes
<b>Optional N2O</b>	No
<b>Emergency O2 Flush</b>	Yes

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights