

BURTONS RAIL MOUNTED ANAESTHESIA MACHINE

102-088-BUNDLE



PRODUCT DESCRIPTION

Fully utilise medical rail systems with this specially designed Burton's anaesthetic machine!

Ideal for piped gas systems this medical rail mounted, O₂ only, anaesthetic machine is a perfect addition to any practice. The machine is built to last, with a stainless steel frame and Burton's own high-quality components. Supplied as standard with a high-accuracy rotameter and emergency O₂ flush (including pressure release safety valve).

Installation and moving the system is easy on your medical rail system, just clip on and go, with built-in lock down points to stop any movement whilst the system is in use.

Use this machine alongside a piped gas system to instantly save precious floor space and remove the need to move large cylinders around the practice.

All models are supplied with :

- 1-metre gas supply lead with indexed BS Schrader probes
- High-accuracy rotameters for precise control of gas flow rates
- Corner unit outlet with built-in emergency oxygen flush
- Non-return valve and safety valve
- 304 stainless steel frame - High-quality, medical grade frame - durable, robust, easy-clean and manufactured to last

Find out more about the key features and benefits of a Burton's Anaesthesia Machine in our Buyer's Guide here:

[Burtons | Choosing an Anaesthesia Machine for your practice](#)

Technical Specification:

- O₂ only
- Overall size: 21cmD x 40cmW x 47cmH

- Outlet connector: BS22mm Male
- Flow control: 0.1 to 10 l/min

Easy Installation

Each unit is supplied with a universal fitting compatible with the standard MediRail system. Our Rail-Ready systems come with universal clamps, allowing the unit to be easily lifted on and off the rail, making it simple to move between locations as needed.

We've purposefully incorporated this design into our product range, recognising that many practices already have MediRails of various lengths installed. By making our units compatible, we help you maximise space and streamline your setup without the need for additional mounting systems.

Don't forget your Vaporiser

The Series 3 is a temperature compensating, out of circuit vaporiser in a traditional design, familiar to most veterinary practitioners. Calibrated specifically for the specified anaesthetic drug using a laser refractometer. *Please note: Vaporiser Sold Separately*

Fitting Type

Choose from the standard taper mounted fitting (cage mount) (UNI) or the quick release (QR) option; allowing you to easily swap the vaporiser to different anaesthetic machines throughout the practice. Additional quick release compatibility blocks will be required for each anaesthesia machine. (Available separately)

Key filled or Funnel (Pour fill)

Key filled (KF) systems require the addition of key filler adaptors but are generally safer and have a lot less wastage compared to standard funnel filled (FF) vaporisers, saving you money in the long run.

Manufactured to last

Burtons Veterinary - At the heart of veterinary anaesthesia for over 30 years.

We don't just supply it! It was designed and built by Burtons in the UK, made to the exacting high standards of British manufacture since 1981, at our principal manufacturing facility in Staplehurst, Kent. Using the latest CNCs and laser manufacturing processes and high-quality stainless-steel materials, this product has been engineered from the ground up to last!

This allows us to reduce our carbon footprint with incredibly low product mileage, made with fully recyclable materials, which contribute to a circular economy.

Free Expert Commissioning with Every Anaesthesia Machine

At Burtons, we support you from day one. That's why we offer free onsite commissioning with every anaesthesia machine purchase.

Once your order is confirmed, one of our expert service engineers will be in touch to arrange a visit at a time that suits you. During the visit, they will:

- Fit and calibrate your vaporiser
- Test and check that your machine is working perfectly
- Install your concentrator if required
- Answer any questions and walk you through the setup and operation

We'll make sure you're fully set up and ready to go, with the option of additional training from your local product expert if needed.

Please note: bracket fitting or any construction work is not included as part of commissioning.

Annual Anaesthesia Machine Service Plans Available

Keep your machine performing at its best with our annual service plans, designed to maintain accuracy and reliability at the heart of your anaesthesia system. Enjoy peace of mind with ongoing support from our UK-based expert service team, just a phone call away.

Discover your perfect Anaesthesia Machine with Burtons using the 8-page Equipment Guide in the latest edition of the

Burtons Insight Magazine (Pages 12-19) here:

[Burtons Insight \(Edition 6\) - Anaesthesia Machine Equipment Guide](#)

ADDITIONAL INFORMATION

Animal	Mixed
SKU	102-088-bundle
Colour/Style	Stainless Steel
Gas Type	O ²
Gas Supply	'E' Size Cylinder (s), 'F' Size Cylinder(s), Piped Gas
Flowmeter Type	Graduated Rotameter
Concentrator	No
Optional O2 Alarm	Yes
Optional N2O	No
Emergency O2 Flush	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights