

## COMPLETE VETSEAL ELECTROSURGERY UNIT

END-VGMS-BUNDLE



### PRODUCT DESCRIPTION

## The only generator with a fully autoclavable bipolar thermo fusion hand piece

The VetSEAL is an advanced electrocautery generator from the human market, with a fully autoclavable bipolar thermos fusion hand piece for sealing vessels up to 7mm. Vessels are permanently sealed without leaving any foreign material, such as sutures, in the body. A seal is achieved when bipolar current is passed through the tissue, causing the collagen and elastins to fuse without carbonisation. Thermo-fusion of vessels with the VetSEAL is fast and secure and can be repeated as many times as necessary during the procedure.

Complete Kit Includes: VetSEAL unit, single bipolar/vetSEAL footswitch, VetSEAL thermocut cable & hand piece, 5 VetSEAL blades.

### ADDITIONAL INFORMATION

<b>SKU</b>	END-VGMS-bundle
<b>Animal</b>	Mixed
<b>Generic Size</b>	Medium

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights