

HUGER 68HS VIDEO GASTROSCOPE - CLEARANCE

END-HGVG-68HS-CL



PRODUCT DESCRIPTION

HUGER AGVE-68 series GI endoscopes employ advanced CMOS image sensor and combined with the compact VIS-68 video Huger endoscopy systems (with superior algorithm technology and refined functions) present exceptionally high-resolution images with clarity in every detail.

- 1240*720p CMOS high resolution endoscope
- New generation handle, smaller, lighter and of remote control button and enhanced maneuverability
- Smoother, flatter and waterproof endoscope body for easier cleaning and disinfection

S/N: AGHS682192

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-HGVG-68HS-cl
Colour/Style	Black
Brand	Huger
Procedure Type	Bronchoscopy, Cystoscopy, Respiratory Examination, Rhinoscopy

Biopsy Channel	2.8mm
Working Length (mm)	1500
Outer Diameter (mm)	8.0

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights