

HUGER 68HS VIDEO GASTROSCOPE - CLEARANCE

END-HGVG-68HS-CL



PRODUCT DESCRIPTION

HUGER AGVE-68 series GI endoscopes employ advanced CMOS image sensor and combined with the compact VIS-68 video Huger endoscopy systems (with superior algorithm technology and refined functions) present exceptionally high-resolution images with clarity in every detail.

- 1240*720p CMOS high resolution endoscope
- New generation handle, smaller, lighter and of remote control button and enhanced maneuverability
- Smoother, flatter and waterproof endoscope body for easier cleaning and disinfection

S/N: AGHS682192

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-HGVG-68HS-cl
Colour/Style	Black
Brand	Huger
Biopsy Channel	2.8mm

Working Length (mm)

1500

Outer Diameter (mm)

8.0

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights