

## COUPLANDS ELEVATORS

128-COUPLANDS



## PRODUCT DESCRIPTION

- Thin cross section at tip to allow good access & positioning on the tooth root.
- High mechanical advantage (ratio of handle diameter to tip width)
- Large & Shaped "American" style handle for comfortable palm grip
- Stainless Steel surgical construction Used to: loosen a tooth root as part of tooth extraction. Can also be used to cut the epithelial attachment to the tooth.

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	128-couplands
<b>Brand</b>	Big O
<b>Dental Instrument Type</b>	Elevators

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

