

## EXAMION PX40 BT PLUS BATTERY GENERATOR

XR-EX332PX40BT



### PRODUCT DESCRIPTION

---

## High-Frequency X-Ray Generator

Perfect choice for optimum power

- **Max. output:** 2,4 kW
- **Tube:** Toshiba D - 124
- **Output:** 100 kV / 35 mA
- **Weight:** 14.2
- **Adjustable kV:** 40 - 100 kV - (1 kV steps)
- **mAs range:** 0.4 - 100 mAs (35 steps)
- **Current range** mA output: 16 - 35 mA
- 35mA @ 40~60kV/0.4~100mAs
- 30mA @ 61~70kV/0.4~64mAs
- 25mA @ 61~70kV/80~100mAs
- 30mA @ 71~80kV/0.4~32mAs
- 25mA @ 71~80kV/36~80mAs
- 25mA @ 81~90kV/0.4~40mAs
- 20mA @ 81~90kV/45~80mAs
- 20mA @ 91~100kV/0.4~50mAs
- 16mA @ 91~100kV/64~80mAs
- **Focal spot:** 1.2 x 1.2 mm
- **Total filtration:** 2.2 mm Al (including collimator)
- **Collimator light field:** Max: 35 x 35 cm @66 cm SID / 40 x 40 cm @ 75 cm SID Min: 5 x 5 cm @ 100 cm SID
- **Collimator lamp Duration:** 30 sec
- **SID measurement:** dual laser pointer & tape measure
- **Dimensions (L x W x H) mm:** 349 x 250 x 214 mm

### Battery

- **Type:** Li-polymer
- **Capacity (mAh):** 4800
- **Exposures (approx.):** 200 On low exposure volume up to one month per charge
- **Charging time:** Less than 3 hours

#### Included accessories

- Memory exposures ( 21 APR standard)
- Aluminum carrying case with shoulder strap
- Podoblok/ MXR accessory compatible
- Optional: EXAMION® - VET Alu X-ray table Professional with sliding Bucky/Column

#### Warranty

- 12 months on hardware / 24 months on software

#### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	XR-EX332PX40BT
<b>Brand</b>	Examion

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights