

## LFM MAGNIFIED EXAM LIGHT

145-810-BUNDLE



### PRODUCT DESCRIPTION

## Shadow-free magnification

The 127 mm glass lens is positioned in the middle of the circular tube, the ideal way of illuminating an object, virtually eliminating shadows

- 3 dioptre lens gives 175% magnification.
- Spring counterbalanced for stability.
- Hinged dust cover for lens.
- Supplied with table mounting bracket.

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	145-810-bundle

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

