

PIN INDEX MULTI-STAGE MA7 REGULATOR 7 BAR

102-224

**PRODUCT DESCRIPTION****Multi-Stage Pin Index Medical Air Regulator with Schrader**

This high quality Burtons multi stage regulator is suitable for all pin index connections to 'E' and 'J' size cylinders with 0-10 bar outlet pressure.

- Pin Index
- Suitable for 'E' and 'J' cylinders
- 0-10 bar pressure

ADDITIONAL INFORMATION

Animal	Mixed
SKU	102-224
Cylinder Size	J, E
Gas Type	MA7

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights