


## PODOBLOCKS

XR-EXTB-PODO



## PRODUCT DESCRIPTION

Select Your Topblock with your order

### Key Features:

- Made of aluminium
- Weight Podoblock: 6,5 kg
- 4 rubber feet
- 120mm
- LM / AP / 45° Handsfree system
- To be used in combination with our Podoholder, Podo-Protectionbox or Cowboybox, for optimal radiation distance
- Avoids fingers on your X-ray

## ADDITIONAL INFORMATION

<b>SKU</b>	XR-EXTB-PODO
<b>Animal</b>	Mixed
<b>Colour/Style</b>	Black

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights