

EXAMION X-R STATIC BT & WS

XR-EXBTWSV



PRODUCT DESCRIPTION

Includes:

- BT 4 way Patient Table
- X-R More Static WS Bucky Wall Stand
- X-R More Static Floor Mounted Tube Stand
- X-R More Generator 50kW CP
- X-Ray Generator CPI CMP200DR
- X-HW More Handswitch, for prep & Exposure
- X-R Tube E7252X TO
- X-R More Stator Cable 12m
- 2 X-R More High Voltage Cable 10m
- X-R More Collimator RO

Bucky table with floor mounted tube stand

Table height - 692 mm (27.2")

Tabletop dimension - Standard: 2260 x 755 mm (89" x 30")
Optional: 2000 x 755 mm (79" x 30")

Tabletop mobility - Longitudinal movement: ± 470 mm (18.5");
optional ± 340 mm at 2000 mm (13.4" at 79")
Transverse movement: ± 150 mm (6")

Object - film/detector - distance - 67 mm (2.6")

Al equivalent value - @ 100 kV: 0.9 mm (0.04") Al; optional 0.7 mm (0.03") Al

Patient bearing load - max.: 230 kg (507 lbs); optional 250 kg (551 lbs)

Focal spot vertical travel - horizontal X-ray beam - 250 - 1892 mm (9.8" - 74.5")

Focal spot vertical - film distance - max. 1267 mm (49.9")

Focal spot vertical - tabletop distance - max. 1200 mm (47.2")

Angulation X-ray tube assembly - Around horizontal support arm: $\pm 180^\circ$
Around tube column (optional): $\pm 90^\circ$
Detents at: $- 90^\circ, 0^\circ, + 90^\circ$

Vertical travel supporting arm - 1642 mm (64.6")

Tube stand longitudinal travel - 1285 mm (50.6")

Bucky wall stand

Vertical movement of the acquisition unit - 370 - 1680 mm (14.6" x 66.1")

Column height - 2133 mm (84")

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	XR-EXBTWSV
Brand	Examion
System Type	Floor Standing

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights