

BULLNOSE OXYGEN SINGLE STAGE REGULATOR 4 BAR WITH SCHRADER

102-230



PRODUCT DESCRIPTION

Bullnose oxygen regulator with schrader

Bullnose oxygen single stage regulator with schrader valve outlet and 4 bar outlet pressure for 'F' and 'G' sized cylinders.

- Oxygen
- Schrader valve
- 4 bar outlet
- Fits cylinders 'F' and 'G' sized

ADDITIONAL INFORMATION

Animal	Mixed
SKU	102-230
Cylinder Size	F, G
Gas Type	O ²

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights