

BURTONS VET FLUID WARMER

129-310



PRODUCT DESCRIPTION

Our new and improved fluid warmer is compact and light weight and now comes with a colourful LCD screen for complete ease of use

Key Features:

- Compact and light weight - 550g
- Colourful LCD screen with detailed info
- Double groove heating for two infusion sets, suitable for fluid and blood
- Fast heating, 3~5 minutes to reach the set temperature
- Cut off power supply automatically if fluid is overheated
- With over-temperature & low temperature alarm
- Two flexible installations: with hanger or pole clamp
- Compatible with general infusion set/tube, no additional disposables
- Optional working with drop sensor - it will give alarm if there's no drops, suitable with gravity infusion

Specifications

- Tube Size: 3~4.5mm
- Heating Range: 30 °C ~ 50 °C, with increment of 1°C or 0.1 °C
- Accuracy: $\leq \pm 0.5^\circ\text{C}$
- Classification: Class I / type B
- Power: AC 100-240V, 50Hz/60Hz
- Net Weight: 550g
- Dimensions: 172*92*55 mm

ADDITIONAL INFORMATION

Animal	Mixed
SKU	129-310
Brand	Burtons
Screen Type	Digital Display
Flow Rate	100-300ml/h

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights