

BURTONS FOLD DOWN CONSULT WALL MOUNTED TABLE

225-215



PRODUCT DESCRIPTION

This economy fold down table is a simple solution for installing a Consult table into a confined area

Its high quality stainless steel frame offers a stable platform for many patients up to 80kg and when your finished it is easy and quick to fold away. Constructed from a fully welded stainless steel frame that is easily cleaned and mounts the table securely to the wall allowing the table to hold up to xxkg*.

- Table top size: 610mm x 885mm
- Complete with black mat
- Maximum load; 80Kg, *when evenly distributed and correctly mounted to wall.

Other sizes available to special order. Please call for a quote

Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	225-215

Colour/Style

Stainless Steel

Brand

Burtons

TiltingTop

No

Mobile

No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights