

BURTONS CIRCLE CIRCUIT SEMI-DISPOSABLE

103-086



PRODUCT DESCRIPTION

A simple circle circuit suitable for patient weights 10kg and above

With both inlet and outlet ports, save gas by recycling expired gases. The sodalime filled circuit will absorb the CO₂, filtering the gas instead of wasting it.

Suitable for use with out of circuit vaporisers, the canister fits directly onto the outlet of the anaesthetic machine (BS22mm), an extension piece is available as optional extra to enable use on non standard anaesthetic machines.

Manufactured by Burtons In the UK - Bringing you the best quality, at the best price!

Key Features

- MRI compatible
- Refillable cannister
- Holds 400g Sodalime
- Fits most anaesthesia machines with British Standard 22mm connector
- Replaceable bag and patient tubing
- 1.6mtr pallelle tubing set
- Patient range 10-40kg
- Saves gas due to recycling via soda lime

ADDITIONAL INFORMATION

Animal

Small Animal

SKU

103-086

Colour/Style

Clear

MRI Compatible

No

Circuit Type

Circle System

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights