

## NON DISPOSABLE CO-AXIAL BAIN CIRCUIT

103-055



### PRODUCT DESCRIPTION

**A non rebreathing circuit with co-axial tubing; fresh gas is supplied through the centre bore whilst exhaust gases are exhaled through the outer tubing**

Recommended for use with patients above 10Kg body weight.

- Suggested minimum flow rate 2.5 times minute volume
- 1.6 metre Co-axial tubing set
- 2M/15F patient connector
- 30mm scavenge valve
- 2 litre rebreathing bag
- Typically for patients > 10Kg
- Safety factor 2.5 x minute volume
- Mapleson D Classification

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	103-055
<b>Durability</b>	Non-Disposable
<b>Tubing Length</b>	1.6 Meter
<b>MRI Compatible</b>	No
<b>Circuit Type</b>	Bain

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights