

V SOCKET MULTI MOUNTING BLOCK KIT

103-145



PRODUCT DESCRIPTION

V Socket Multi Mounting Block Kit

The black V Socket Multi Mounting Block Kit enables you to mount an Air Break or Cyclo-Flo directly to a Burtons anaesthetic machine, clamping onto the box section of any Burtons anaesthetic machine. This is via the two threads, that align with the corresponding back plate, bolting into place.

The socket can also be used for mounting directly onto a suitable wall, using suitable fixings for the desired wall (not supplied) that can hold the weight of the unit, using the fixing points that sit centrally within the recess of the 'V' block making it easy for you to attach to the wall.

Suitable for securing Burtons Cyclo-Flo, Cyclo-Flo V2 and Air Breaks.

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	103-145

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights