

OSTER CRYOGENX BLADE, SIZE 40

OST-BLADE-40

**PRODUCT DESCRIPTION**

Oster® Professional blades are produced from the finest high-carbon steel and handcrafted to our rigorous standards, providing unparalleled precision and performance

Oster blades are hardened to RC62-66 (Rockwell Scale C) in-house prior to receiving a cryogenic treatment which further enhances the hardening process. The blade cutting surfaces are then sharpened, finished and tested with proven Oster® proprietary processes, ensuring ultra-sharp blades that continue to stand the test of time.

Key Features:

- 0.25mm

ADDITIONAL INFORMATION

SKU	OST-BLADE-40
Colour/Style	Silver
Blade Size	40
Blade Coating	N/A
Blades Fits	Oster
Blade Type	CryogenX
Suitable for cat or dog	Dog

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights