

## SET OF 24 RED RUBBER ENDOTRACHEAL TUBES

104-060



### PRODUCT DESCRIPTION

## Quality sterilisable Magill endotracheal tubes complete with balloon

#### Key features:

- Cuffed and uncuffed option
- Plain design
- Magill style
- Autoclavable to 134°C
- Reusable
- Do not dry heat

#### Set Includes:

- 1 each of the cuffed size 3.0,3.5,4.0,4.5,5.0,5.5,6.0,6.5,7.0,7.5,8.0,8.5,9.0,9.5,10.0,11.0
- 2 each of uncuffed size 2.5,3.0,3.5,4.0

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	104-060
<b>E.T Tube Size</b>	Full Set
<b>Tube Material</b>	Rubber

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights