

**TOFFELN AKTIVFLEX TRAINER**

220-CONFIG\*

**PRODUCT DESCRIPTION**

**The UK-first EVA enclosed trainer, designed for the healthcare environment to deliver all of the wear with none of the tear.**

A super durable, superior lightweight nursing shoe complete with Grip-Safe® sole and removable insole. Ergonomically designed with input from the University of Salford to deliver impressive performance characteristics.

- Elasticated gusset
- Slip on with ease and enjoy a secure fit
- ComfortTech Insole Included
- Walk with greater comfort with the specially designed insole
- Grip-Safe® Sole
- This unique sole gives excellent slip resistance, and a much greater penetration resistance.
- Side Ventilation
- Allows your feet to breathe, keeping them cool throughout the day
- Stylish Rubber Laces

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	220-config*
<b>Category Sub Heading</b>	Available in 4 colours & all sizes

Male/Female

Unisex

Brand

Toffeln

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights