

**SET OF 6 VETERINARY CLEAR PERSPEX MASKS WITH DIAPHRAGM**

104-223

**PRODUCT DESCRIPTION****Polycarbonate high visibility masks with soft silicone rubber diaphragm to ensure a proper seal**

Set of 6 veterinary mask clear perspex complete with diaphragm. largest sizes are supplied with adjustable harness

**Set Includes:**

X-Large ( 5.1/4in Dia X 5in Deep ) - adjustable harness

Large ( 4 1/4in Dia X 4in Deep ) - adjustable harness

Medium ( 3 1/2in DIA X 3 1/4in Deep ) - adjustable harness

Small ( 2 1/2inDia X 2in Deep ) - adjustable harness

Extra-Small ( 2 1/4in Dia X 1 1/4in Deep )

Tiny ( 1 3/4in Dia X 1in Deep )

**ADDITIONAL INFORMATION****Animal**

Small Animal

**SKU**

104-223

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights