

SET OF 6 VETERINARY CLEAR PERSPEX MASKS WITH DIAPHRAGM

104-223

**PRODUCT DESCRIPTION****Polycarbonate high visibility masks with soft silicone rubber diaphragm to ensure a proper seal**

Set of 6 veterinary mask clear perspex complete with diaphragm. largest sizes are supplied with adjustable harness

Set Includes:

X-Large (5.1/4in Dia X 5in Deep) - adjustable harness

Large (4 1/4in Dia X 4in Deep) - adjustable harness

Medium (3 1/2in DIA X 3 1/4in Deep) - adjustable harness

Small (2 1/2inDia X 2in Deep) - adjustable harness

Extra-Small (2 1/4in Dia X 1 1/4in Deep)

Tiny (1 3/4in Dia X 1in Deep)

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	104-223

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights