

**LITTMANN LIGHTWEIGHT NURSE STETHOSCOPE**

117-671-CONFIG1\*

**PRODUCT DESCRIPTION**

**The 3M™ Littmann® Classic II S.E. Stethoscope is a dependable, hard working clinical tool for physical assessment and monitoring of your patients. It's the choice of millions of care providers who are dedicated to achieving the best patient outcomes.**

**Key Features**

- Versatile chestpiece design provides high acoustic sensitivity with both a tunable diaphragm and open bell.
- Tunable diaphragm: Hear high or low frequency sounds by slightly adjusting pressure on the chestpiece.
- Non-chill rim and diaphragm provide patient comfort.
- Soft-sealing eartips provide an excellent acoustic seal and comfortable fit.
- Headset is easily adjusted for individual fit and comfort. Angled eartubes align with ear canals.

**ADDITIONAL INFORMATION****Animal**

Small Animal

**SKU**

117-671-config1\*

**Stethoscope Type**

Lightweight SE

**Chestpiece Style**

Standard

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights