

## LITTMANN CLASSIC III STETHOSCOPE - BLACK CHESTPIECE

117-670-BK-CONFIG\*



### PRODUCT DESCRIPTION

**The Littmann Classic III stethoscope offers high acoustic sensitivity for exceptional performance. It features dual turn-able diaphragms for adult and paediatric use**

#### Features

- Next-generation tubing provides longer life due to its improved resistance to skin oils and alcohol, making it easier to clean.
- The single-piece turn able diaphragm is easy to attach, and easier to clean due to the smooth surface.
- Paediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.
- Small turn able diaphragm is useful for pediatric, small, or thin patients; around bandages; and for carotid assessment.
- Stainless steel chest piece is precision-machined into an aesthetically pleasing, less angular shape.
- Snap-tight ear-tips have a soft smooth surface for a comfortable fit.
- Five Year Warranty.

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	117-670-BK-config*
<b>Catagory Sub Heading</b>	Colour Options
<b>Stethoscope Type</b>	Classic III
<b>Chestpiece Style</b>	Black

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights