

LITTMANN CLASSIC III STETHOSCOPE - TURQUOISE TUBING - PINK STEM - SMOKE CHESTPIECE

117-670-T-S-P



PRODUCT DESCRIPTION

The Littmann Classic III stethoscope offers high acoustic sensitivity for exceptional performance. It features dual turn-able diaphragms for adult and paediatric use

Features

- Next-generation tubing provides longer life due to its improved resistance to skin oils and alcohol, making it easier to clean.
- The single-piece turn able diaphragm is easy to attach, and easier to clean due to the smooth surface.
- Paediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.
- Small turn able diaphragm is useful for pediatric, small, or thin patients; around bandages; and for carotid assessment.
- Stainless steel chest piece is precision-machined into an aesthetically pleasing, less angular shape.
- Snap-tight ear-tips have a soft smooth surface for a comfortable fit.
- Five Year Warranty.

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	117-670-T-S-P
Colour/Style	Turquoise
Catagory Sub Heading	Colour Options
Stethoscope Type	Classic III
Chestpiece Style	Limited Edition Smoke

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights