

## BURTONS VETAIR CLASSIC MOBILE DENTAL UNIT

118-704-2



### PRODUCT DESCRIPTION

## A compact, all-encompassing mobile dental unit

The Burton's VETair Classic's array of features has been designed with vets and vet nurses in mind, making dental procedures as efficient and straightforward as possible. Whilst its minimal maintenance requirements make it an extremely cost-effective dental system.

Designed & manufactured by Burton's in the UK, the VETair Classic is built upon a 'Lifetime' stainless steel frame, designed to provide ultimate longevity and stability, whilst still being easy to clean. Its compact size means you have everything to hand without taking up extra space that could be better utilised. An accessible storage drawer enables you to keep all your instruments and spares within the dental system for easy access during procedures.

Strong lockable castors enable easy moveability for getting as close to your patient as possible, or for moving it out of the way when not in use. All Burton's dental machines are designed alongside Burton's knee-space tub tables, allowing you to easily tuck your machine away underneath the table when not in use.

Our latest improvement to the VETair Classic is the oil-free compressor, which eliminates the risk of oil contaminants within the dental machine and reduces maintenance costs and downtime. But unlike other oil-free compressors, the VETair removes the stigma of oil-free = noisy. Using Burton's very own noise-dampening technology, which you may have already experienced within our concentrators, this combination of spring-mounted components softens the inertia of starting and stopping the pump, providing a gentle, quiet compressor action. Reducing stress for practice staff.

**Additional starter kit and scaler - adding these will mean you have everything you need to get started, including:** A high-speed LED handpiece for drilling, a low-speed LED handpiece for polishing and finishing and an air/water spray for cleaning away debris and cooling when high-speed instruments are in use.

## Key Features

- Burtons' noise-dampening technology removes the stigma of oil-free meaning noisy.
- Accessible drawer for handpiece and accessory storage
- High-quality, lockable castors for easy movability
- Three station control unit including high speed, low speed & air/water mist.
- The VETair works via a foot-pedal and simple flick switch between handpieces, freeing up hands during procedures.
- A good-sized onboard water bottle means there is no need to be stopping to re-fill water bottles during procedures. This one bottle is used for an additional scaler as well as handpieces for less maintenance of water reservoirs.
- Reduce treatment time to reduce anaesthesia risk by having everything you need in one place.
- Fits underneath Burtons' knee space tub tables when not in use.
- Coiled cables keep trailing leads from getting run over on the floor and help keep the system tidy when not in use.
- A fixed power bank to the back of the machine gives you the flexibility of plugging in additional accessories needed, this is fixed high on the machine, out of the way of any liquid spillages.
- Designed & manufactured in the UK.
- Fully supported by Burtons after-care support team.
- Consumables such as 'O' rings and water bottles can be easily purchased and replaced if, and when, needed.

## Technical Specs:

- Dimensions: 750mmH x 670mmW x 450mmD (including handles and handpiece holders)
- Weight: 52kg
- Power: 240v AC, 50/60hz, 3.6 Amp, 500W
- Volume: 55 dB(A)
- Reservoir capacity: 10ltr)
- Water bottle capacity: 2L

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	118-704-2
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights