

E3 EMS COMPATIBLE SCALER - CLEARANCE

149-000-CL



PRODUCT DESCRIPTION

Economy EMS Compatible Dental Scaler is a high-quality, low-cost dental scaler that offers all the features and benefits of a more expensive unit.

Key Features:

- The handpiece is made of titanium and peek material. Lightweight and ergonomic, autoclavable up to 134°C
- The intelligent tracking and power feedback system provides maximum protection for the gums.
- Multi-functions of scaling, perio and endo, to meet different demands of ultrasonic scaler treatment.
- Touch panel: Simplified operation and precise setting.
- Tips & handpiece compatible with EMS series
- Supplied with 6 tips

Model	E3 Scaler
Mode:	Continuous operation
External voltage	220-240 VAC
Mains frequency	50/60 Hz
Water supply	0.1 to 0.5 bar (100 to 500 kPa)
Frequency range:	24-33 kHz

As with the Piezon scalers, the E3 Scaler can also be fitted with the uniquely developed urinary tips which are inserted into the

urethra to pulverise wedged calculus, a particularly common and recurring problem with domestic cats. The procedure takes only a few seconds and eliminates the potential need to carry out a urethrotomy. Please see the related products for more information.

Note: If this product is not used in conjunction with the VETair dental system you will require a separate water bottle part number 118-015

Clearance products may be non-pristine, for example:

- On occasion some items may show signs of minimal use, which may include minor cosmetic marks which do not affect the functionality in any way.
- Although most of our clearance products are sold in their original manufacturers packaging, there may be occasions, due to packaging damage, where plain white or brown packaging may be used

Clearance units sold with 6 month standard guarantee, unless otherwise specified. This does not affect your statutory rights.

Please see our terms and conditions for more information

S/N: E3K70019

ADDITIONAL INFORMATION

Animal	Mixed
SKU	149-000-cl
Generic Size	Small
Colour/Style	White
Examination Type	Oral Examination, Small Furrries, Tooth Extraction
Fitting Type	4 Hole

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

