

## BURTONS ADAPTABLE2 ANAESTHESIA MACHINE

102-060-BUNDLE\*



### PRODUCT DESCRIPTION

---

## The same industry-leading Adaptable... in a new, smarter material!

**The Adaptable2 is manufactured in a solid uPVC material frame that is impervious to liquid and resistant to bacteria, making for a completely sealed, easy-wipe and waterproof anaesthesia machine that stands up to rigorous cleaning regimes and reduces contamination risk.**

Built upon the original Burton's AAA (Adaptable Anaesthesia Apparatus) frame concept, the Adaptable2 Anaesthesia System is designed for versatility and can form the basis of a variety of anaesthesia circuits. The Triple A is a clever, convenient and modular frame, designed specifically for the veterinary market, to accommodate bespoke configurations. Perfect for piped gas, this machine can also accommodate two F-size oxygen cylinders in hinged clamps on the rear and an onboard oxygen concentrator on the base shelf. Any or all of these oxygen supply methods can be used without the need to reconfigure your machine's connectors.

For peace of mind in theatre, the Adaptable2 includes a number of carefully considered safety features. In cases of power outage, the manual changeover valve allows you to switch from your onboard concentrator quickly and seamlessly to your backup cylinders. The emergency flush button contains two integral safety valves for pressure relief and to prevent expired gases from flowing back through the circuit to the vaporiser and beyond. Safety is top priority with all Burton's machines having a built-in emergency O<sub>2</sub> flush with pressure safety valve as standard.

### Key Features

- Frame made from a waterproof, antibacterial, solid uPVC – a continuous easy-wipe surface
- High-accuracy 0-10l/min logarithmic rotameter;
- Emergency O<sub>2</sub> flush includes integral pressure relief and non-return safety valves;
- Connects seamlessly with the piped oxygen in your practice, F-size cylinders, or with an on-board concentrator;
- Rapid manual change-over valve for switching between concentrator and cylinder oxygen supplies – vital for power outages;
- Includes a monitor shelf, situated on the top of the frame;
- Can accommodate up to two vaporisers;
- Fully mobile, solid four-wheel base with independent wheel locks;
- Full warranty, supported by Burtons' nationwide Service Team;
- Designed and manufactured by Burtons in the UK
- This machine frame is 1380mm High x 480mm Wide x 650mm Deep. With the optional Cyclo-Flo & Air Brake fitted the width increases to 630mm (See dimensions diagram above)
- F-size cylinder rings and regulator kit, concentrator and vaporiser(s) sold separately.

Find out more about the key features and benefits of a Burtons Anaesthesia Machine in our [Buyer's Guide](#).

#### Technical Specification:

- Dimensions: 1380 x 480 x 650mm HWD (see diagram of dimensions above)
- Weight: 20kg

#### Don't forget your Vaporiser

The Series 3 is a temperature compensating, out of circuit vaporiser in a traditional design, familiar to most veterinary practitioners. Calibrated specifically for the specified anaesthetic drug using a laser refractometer. *Please note: Vaporiser Sold Separately*

#### Fitting Type

Choose from the standard taper mounted fitting (cage mount) (UNI) or the quick release (QR) option; allowing you to easily swap the vaporiser to different anaesthetic machines throughout the practice. Additional quick release compatibility blocks will be required for each anaesthesia machine. (Available separately)

#### Key filled or Funnel (Pour fill)

Key filled (KF) systems require the addition of key filler adaptors but are generally safer and have a lot less wastage compared to standard funnel filled (FF) vaporisers, saving you money in the long run.

#### Manufactured to last

Burtons Veterinary - At the heart of veterinary anaesthesia for over 30 years.

We don't just supply it! It was designed and built by Burtons in the UK, made to the exacting high standards of British manufacture since 1981, at our principal manufacturing facility in Staplehurst, Kent. Using the latest CNCs and laser manufacturing processes and high-quality stainless-steel materials, this product has been engineered from the ground up to last!

This allows us to reduce our carbon footprint with incredibly low product mileage, made with fully recyclable materials, which contribute to a circular economy.

#### Free Expert Commissioning with Every Anaesthesia Machine

At Burtons, we support you from day one. That's why we offer free onsite commissioning with every anaesthesia machine purchase.

Once your order is confirmed, one of our expert service engineers will be in touch to arrange a visit at a time that suits you. During the visit, they will:

- Fit and calibrate your vaporiser
- Test and check that your machine is working perfectly
- Install your concentrator if required
- Answer any questions and walk you through the setup and operation

We'll make sure you're fully set up and ready to go, with the option of additional training from your local product expert if needed.

*Please note: bracket fitting or any construction work is not included as part of commissioning.*

### Annual Anaesthesia Machine Service Plans Available

Keep your machine performing at its best with our annual service plans, designed to maintain accuracy and reliability at the heart of your anaesthesia system. Enjoy peace of mind with ongoing support from our UK-based expert service team, just a phone call away.

**Discover your perfect Anaesthesia Machine with Burtons using the 8-page Equipment Guide in the latest edition of the Burtons Insight Magazine (Pages 12-19) here: [Burtons Insight \(Edition 6\) - Anaesthesia Machine Equipment Guide](#)**

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	102-060-bundle*
<b>Category Sub Heading</b>	Configure to your practice needs!
<b>Gas Type</b>	O <sup>2</sup>
<b>Gas Supply</b>	'F' Size Cylinder(s)
<b>Concentrator</b>	No
<b>Optional O2 Alarm</b>	No
<b>Optional N2O</b>	No
<b>Emergency O2 Flush</b>	No

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights